

## BD LSRFortessa Cell Analyzer (5532)

	<b>Detector</b>	<b>Mirror</b>	<b>Filter</b>	<b>*Suggested Fluors</b>
<b>Blue (488nm)</b>	FITC	505LP	530/30	BB515, FITC, FVS520
	PerCP-Cy5.5	690LP	710/50	PerCP-Cy5.5, 7-AAD
<b>Yellow Green (561nm)</b>	PE		586/15	PE, TdTomato
	PE-CF594	600LP	610/20	FVS620, mCherry
	PE-Cy5	635LP	670/30	
	PE-Cy7	750LP	780/60	PE-Cy7
<b>RED (640nm)</b>	APC		670/30	APC, Alexa647, Sytox Red
	Alexa 700	690LP	730/45	Alexa700, FVS700
	APC-Cy7	750LP	780/60	Fixable near IR , FVS780
<b>Violet (405nm)</b>	BV421		450/50	BV421, Pacific Blue, V450
	BV510	495LP	525/50	BV510
	BV605	595LP	610/20	BV605
	BV650	635LP	670/30	BV650
	BV711	690LP	710/50	BV711
	BV786	750LP	780/60	BV786
<b>UV (355nm)</b>	BUV396		379/28	FVS440UV, BUV396
	BUV496	450LP	515/30	Fixable Aqua, DAPI
	BUV805	770LP	820/60	BUV805